

StaebeL Deluxe Foam Cushion User Manual

Product Code: 321BC



Thank you for choosing Staelbel Deluxe Foam Cushion

We aim to ensure you are satisfied with your cushion in terms of both performance and durability. We therefore kindly ask you to spend a few minutes reading through these user instructions to ensure you can take the best possible care of your cushion. If you have any questions that are not answered in these user instructions, please contact us at Staelbel. You will find our contact details on the user instructions.

- The product must only be used for its intended purpose.
- DO NOT USE in the vicinity of or in contact with fire sources/hot surfaces, such as fire, burning cigarettes or hot lamps.
- Position and use the product in a manner that prevents it from becoming trapped or damaged.
- Always make sure the cushion is the right size for the patient.
- The cover does not allow liquid or air to penetrate but is vapour permeable.
- Prevent damage to the hygiene cover from sharp objects.

INTRODUCTION

A comfortable cushion for the low risk user, the deluxe pressure relief foam cushion flows under pressure to provide a fluid like support.

- This product has been tested to a safe working load of 100kg.

Specifications

Product Code	3218	3220	3221	3223	3224	3225	3226	3227
Width	400 mm	450 mm	500 mm	460 mm	400 mm	450 mm	500 mm	460 mm
Length	400 mm	400 mm	450 mm	560 mm	510 mm	450 mm	500 mm	610 mm
Thickness	90 mm	90 mm	90 mm	90 mm	90 mm	90 mm	90 mm	90 mm
SWL	100 kg	100 kg	100 kg	100 kg	100 kg	100 kg	100 kg	100 kg

CLEANING/DISINFECTION

Foam Core

- The hygiene cover protects the foam core completely, which means cleaning seldom is necessary. If, for any reason, the core needs to be cleaned, use a mild cleaning solution and water (for example spray on). Press out the water carefully.
- Let the foam core dry in a warm, ventilated area, without direct exposure of sunlight. The foam core needs to be completely dry before it is used again.

Cover

- Some surface wrinkling may result from cleaning procedures. This has no adverse effect on the cushion cover properties
- Abrasive cleaning agents should not be used.
- The cover can be tumble dried provided the operator ensures it is removed as soon as all moisture is removed. The cover must not be left sitting on the heated bowl of the dryer when drying is complete.
- Do not iron the cover
- If customers have their own particular cleaning methods that must be used and that are not covered in this article, they should consult with Staelbel for further direction.

- The cushion Cover needs to be treated with care when cleaning the surface. Many industrial cleaning agents are extremely harsh on the coating and can cause it to break down. Please refer to the list of cleaning products that have been tested and approved for use on the Cover.
- To avoid premature degradation of the surface coating never use an abrasive cloth, sponge or other similar device on the surface. Carefully follow manufacturer instructions for recommended exposure time of cleaner on surface. Wipe off cleaning product residue with a clean, soft cloth after cleaning and allow the surface to dry completely before use.
- The normal cleaning instructions are to wipe the surface clean with a sponge or cloth with warm soapy water.
- The Cover can be washed in warm soapy water at up to 95°C.
- Do not use solvent on the cover.
- The cushion cover is dry-cleanable with all dry cleaning fluids except trichloroethylene.
- The cover is able to be autoclaved at a temperature of 134°C for 3 minutes. However, care must be taken to first clean the surface with pure water and then to make sure cushions are not stacked on top of each other without any air gap separation. Spacers should be placed between cushions when a number of them are being autoclaved at the same time.

STORAGE

Store in a cool dry area. Avoid excessive pressure and contact with non-colourfast materials. Make sure the cushion is stored where it will not be exposed to thrusts, hits or other things that can damage it.

SERVICE & WARRANTY

We warrant the product to be free from manufacturing defects for the period of **12 months** from date of purchase. We will make available upon request; component part lists, descriptions, or other information that will assist service personnel to repair those parts of equipment that are designated by the manufacturer as repairable by service personnel.