


StaebeL Custom Sling Order Form

Thank you for your interest in our StaebeL Patient Slings. Please complete the below steps below by ticking the box for your preference.



Form Information

Full Name: _____ Contact: _____ Date: _____

General Sling Design

Base Sling Design (loop fixing only)

- General Purpose
 General Purpose Head Support
 Care
 Care Head Support
 Junior
 Amputee
 Walking / Standing

Base Sling Material

- Mesh
 Fabric
 Parasilk

Base Sling Type

- Pivot
 Pivot Care
 Standing Transfer
 Super Soft
 Other

	Loop	Clip	E Clip	Dual (Clip)	Dual (E Clip)
<input type="radio"/> Pivot					
<input type="radio"/> Pivot Care					
<input type="radio"/> Standing Transfer					
<input type="radio"/> Super Soft					
<input type="radio"/> Other					

Base Sling Size

- XS (Extra Small)
 S (Small)
 M (Medium)
 L (Large)
 XL (Extra Large)
 XXL (Extra Extra Large)
 Other*

* If you have selected 'other', please provide generic dimensions (A - D) plus any other relevant dimensions in the sketch box provided on the last page.


Sling Modifications

MOD 001 - Change / Add Torso Belt

- Buckle
 Neoprene/Buckle Over
 Custom Length
 Velcro
 Neoprene/ Velcro
 Custom Width

Comments


MOD 002 - Add / Reduce Leg Length



Add Leg Length _____ mm

Reduce Leg Length _____ mm


MOD 003 - Add / Reduce Leg Width




Add Leg Width _____ mm

Reduce Leg Width _____ mm


MOD 004 - Add / Reduce Head Support



MOD 005 - Add / Reduce Body Length



MOD 006 - Add / Reduce Body Width



NB - Body width adjustment will not affect generic external profile of sling.

Add Head Support _____ mm
 Add Body Length _____ mm
 Add Body Width _____ mm
Reduce Head Support _____ mm
 Reduce Body Length _____ mm
 Reduce Body Width _____ mm

MOD 007 - Add Modesty Straps

Modesty Straps keeps legs of sling together without crossing




MOD 008 - Care Sling Leg Strap

Assists in support of upper thighs for clients with low tone using care slings.



MOD 009 - Stand Up Hemi Band

Supports arms affected by hemiplegia in a standing transfer



MOD 010 - Stand Up Drop Strap


Assists clients with limited weight bearing in a standing transfer.



Sling Sketch Box & Dimensions

* If you have selected 'other' in base sling size, please provide generic dimensions (A - D) plus any other relevant dimensions in the sketch box provided below.

Additional Comments:



Sling Sketch and Dimensions

A=

B=

C=

D=

Disclaimer: While Custom Slings are covered by manufacturer's warranty, no responsibility is taken for their suitability and functionality.